

Lab - Measurement & Calculations

Purpose:

Practice making measurements and use dimensional analysis to calculate with appropriate significant figures.

Inquiry:

What is the average thickness of aluminum foil, both heavy duty and regular?

Thickness of a thin sheet is ~~part~~ one dimension of both volume and density.

As Aluminum foil is quite thin, a ruler is not precise enough to determine thickness. I will need to use Density with mass and length & width to calculate the thickness.

Safety:

All general lab safety precautions apply. Aluminum and water are not corrosive or dangerous. Take care when using glass graduated cylinders.

Procedure:

$$(D) \text{ Density} = \frac{\text{mass}(M)}{\text{Volume}(V)} \quad (V) \text{ Volume} = (L) \text{ Length} \times (W) \text{ Width} \times (T) \text{ Thickness}$$

$$D = M \times \frac{1}{L \times W \times T}$$

$$T = \frac{V}{L \times W}$$

$$D = \frac{M}{L \times W \times T}$$

If I can determine the mass, the length the width and the density, I can calculate the thickness.

$$T = \frac{M}{L \times W \times D}$$

Steps

- 1) Use granulated aluminum to determine density by taking a small sample & using a scale to determine mass.
- 2) Use water displacement of granulated aluminum sample to determine volume.
- 3) Obtain multiple trials worth of data. Take an average from each set of measured values.

$$D = \frac{M}{V}$$

Observations

- 1) scale was limited to ± 0.1 g
Aluminum granules appeared slightly oxidized
- 2) graduated cylinder measured to ± 0.1 mL
Some Aluminum dust floated on top.
- 3) third trial yielded a displacement of less than 1 mL - limiting sig figs to 1 place.

Steps

- 4) Obtain two samples each of regular & heavy duty Aluminum foil. Cut foil into a rectangle with 90° corners. Length of both sides should be consistent.
- 5) Measure the length, width, and mass of each sample.
- 6) Use the average density, the length, width, and the mass to calculate thickness for both regular samples & both heavy duty samples.

Data:

Density Determination

Trial 1 (granules)	Mass(g)	3.65	Volume(mL)	1.4
Trial 2	(g)	2.77	(mL)	1.0
Trial 3	(g)	1.99	(mL)	.8

Density

Trial 1	Trial 2	Trial 3
$\frac{3.65g}{1.4cm^3} = 2.6g/cm^3$	$\frac{2.77g}{1.0mL} = 2.8g/cm^3$	$\frac{1.99g}{.8} = 2.5g/cm^3$
Density Average = $\frac{2.6g}{cm^3}$		

Calculations

Regular Foil

Trial 1

$$T = .128g \times \frac{1}{7.57cm} \times \frac{1}{3.45cm} \times \frac{1cm^3}{2.6g} = .001767cm \approx \boxed{1.7 \times 10^{-3}cm}$$

Trial 2

$$T = .15g \left(\frac{1}{7.62cm} \right) \left(\frac{1}{4.15cm} \right) \left(\frac{1cm^3}{2.6g} \right) = .001824cm \approx \boxed{1.8 \times 10^{-3}cm}$$

Average

$$\boxed{1.8 \times 10^{-3}cm}$$

Observations

- 4) Regular & Heavy Duty foil was cut to 90° angles. A pen was used to measure and cut angles - some ink was present.

- 5) Ruler measured to $\pm .01cm$

- 6) Density calculation was limited to 2 sig figs.
Length measurements were limited to 3 sig figs
Mass measurements limited to 2 sig figs.

Length, Width, & Mass

Regular

$$L(cm) \underline{7.57} \quad W(cm) \underline{3.45} \quad (g) \underline{.12}$$

$$L(cm) \underline{7.62} \quad W(cm) \underline{4.15} \quad (g) \underline{.15}$$

Heavy Duty

$$L(cm) \underline{4.96} \quad W(cm) \underline{3.41} \quad (g) \underline{.13}$$

$$\underline{4.78} \quad \underline{4.70} \quad \underline{.16}$$

$$T = M \times \frac{1}{L} \times \frac{1}{W} \times \frac{1}{D}$$

Heavy Duty Foil

Trial 1

$$T = .13g \left(\frac{1}{4.96 \text{ cm}} \right) \left(\frac{1}{3.61 \text{ cm}} \right) \left(\frac{1 \text{ cm}^3}{2.6g} \right) = .002792 \approx 2.8 \times 10^{-3} \text{ cm}$$

Trial 2

$$T = .16g \left(\frac{1}{4.78 \text{ cm}} \right) \left(\frac{1}{4.70 \text{ cm}} \right) \left(\frac{1 \text{ cm}^3}{2.6g} \right) = .002739 \approx 2.7 \times 10^{-3} \text{ cm}$$

Average thickness

$$\boxed{2.8 \times 10^{-3} \text{ cm}}$$

Claim:

- The average thickness of regular foil is $1.8 \times 10^{-3} \text{ cm}$
- The average thickness of heavy duty foil is $2.8 \times 10^{-3} \text{ cm}$

Evidence and Analysis:

Thickness could have been measured by taking a volume measurement and dividing by an area measurement. The extremely thin nature of the foil made volume measurements with available equipment ineffective.

As aluminum is an element and density is an intensive property of aluminum, use of granulated aluminum for volume and mass measurements would transfer to the apply to the foil sample. Calculating the density of aluminum would allow the use of the foil mass in dimensionally determine thickness. (See Detailed Procedure & Calculations)

Reading and Reflection.

According to Elmhurst College, Virtual Chembook (Ophardt, 2003) the accepted value for the density of aluminum is 2.7 g/cm^3 . The calculations we used determined a density of 2.6 g/cm^3 . According to U.S. Packaging and Wrapping, LLC (2019) standard Duty Aluminum foil is measured ~~between .0004 mil and .0007 mil~~ thick. The unit "mil" is $1/1000^{\text{th}}$ of an inch. Standard ~~mil~~ foil is measured between

$$\begin{aligned} .0004 \text{ mil} & \left(\frac{1 \text{ in}}{1000 \text{ mil}} \right) \left(\frac{2.54 \text{ cm}}{1 \text{ in}} \right) = .0004 \text{ m and } .0007 \text{ in} \\ .0004 \text{ in} & \left(\frac{2.54 \text{ cm}}{1 \text{ in}} \right) = 1.0 \times 10^{-3} \text{ cm} \\ .0007 \text{ in} & \left(\frac{2.54 \text{ cm}}{1 \text{ in}} \right) = 1.8 \times 10^{-3} \text{ cm} \end{aligned}$$

the reported values from the packaging company correspond with the values measured in the lab

Heavy Duty Foil is recorded as .001 in and .0008 in thick (U.S. Packaging, 2013)

$$.0008 \text{ in} \left(\frac{2.54 \text{ cm}}{1 \text{ in}} \right) = .0020 \text{ cm} \approx 2.0 \times 10^{-3} \text{ cm}$$

$$.001 \text{ in} \left(\frac{2.54 \text{ cm}}{1 \text{ in}} \right) = 2.5 \times 10^{-3} \text{ cm}$$

The measured lab values for Heavy Duty Foil in the lab slightly exceed reported values.

The precision with which our measured values were reported (limited by the precision of the ruler, scale, and graduated cylinder) tend to correlate with the significant figures observed from both the academic source (Zumdahl, 2013) and the commercial source (U.S. Packaging, 2013).


We may be able to perform the analysis of thickness of aluminum foil with one additional significant figure by increasing the aluminum granule sample to insure water displacement of greater than one mL. However, the data reported as taken, truly limits precision to only 1 sig. fig.

Percent by Mass of NaHCO_3 and KCl Mixture.

Purpose/Questions:

Sodium bicarbonate and potassium chloride are both white powders. When mixed in an unknown proportion, chemical analysis can be used in determination of percent mass of each compound in the mixture.

Sodium bicarbonate will decompose when heated into water vapor and carbon dioxide.



Potassium chloride will not decompose when heated. The difference in boiling points between NaHCO_3 , Na_2CO_3 and KCl are sufficiently high to allow careful heating which will result in decomposition of NaHCO_3 , without affecting melting of any other substance present.

After decomposition of NaHCO_3 , the ratio of CO_2 & H_2O lost to the original amount of NaHCO_3 will allow us to determine initial amount of NaHCO_3 present. Using stoichiometric ratios it is possible to determine the amount of CO_2 & H_2O that will be lost per gram of NaHCO_3 .

$$1.000 \text{ g NaHCO}_3 \left(\frac{1 \text{ mol}}{84.007 \text{ g}} \right) \left(\frac{1 \text{ mol H}_2\text{O}}{2 \text{ mol NaHCO}_3} \right) \left(\frac{18.013 \text{ g}}{1 \text{ mol H}_2\text{O}} \right) = \frac{.107 \text{ g H}_2\text{O}}{\text{g NaHCO}_3}$$

$$1.000 \text{ g NaHCO}_3 \left(\frac{1 \text{ mol}}{84.007 \text{ g}} \right) \left(\frac{1 \text{ mol CO}_2}{2 \text{ mol NaHCO}_3} \right) \left(\frac{44.009 \text{ g}}{1 \text{ mol CO}_2} \right) = .262 \text{ g CO}_2$$

Sodium bicarbonate will decompose at the constant ratio of .369 g CO_2 and H_2O , for every 1.000 g of NaHCO_3 .

By observing the mass lost in a $\text{NaHCO}_3/\text{KCl}$ mixture when heated, the above ratios can be used to determine % NaHCO_3 in the mixture.

Safety Considerations:

Glassware will need to be prepared using HCl which is corrosive to skin and must be neutralized before

- Use of Hot glassware will require eye Protection and careful use of tongs to maneuver.

Procedure/Observations:

Procedure

- 1) Clean porcelain dish with 6M HCl, let soak, rinse thoroughly, and heat to constant mass let cool.
- 2) Mass dish and cover and record. Add approx 1g of KCl/NaHCO₃ mixture and record mass of mixture.
- 3) Heat sample gently for 5 min. Ensure sample does not melt. Heat continuously for 15 min.
- 4) Heat sample to constant mass.

Observations

- 1) After cleaning, small amounts of brown contaminant remained on dish.
- 2) Contents for sample were taken from top of mixture, the mixture was hand mixed.
- 3) After 4 min, edges of sample started to melt. Heat was reduced. Sample formed crusty, cracked, appearance.
- 4) Sample measured within .01g of previous measurement after 5 min of additional heating.

Data Table:

Trial 1	Trial 1
mass crucible	100.09 g
mass w/ mixture	101.44 g
mass mixture	1.35 g
mass after 15 min heat	101.01 g
mass after 5 min cool	101.00 g
mass lost CO ₂ /H ₂ O	.44 g

Trial 2
100.10 g
101.64 g
1.44 g
101.31 g
101.30 g
.54 g

Calculations mass % NaHCO₃:

Trial 1

$$.44 \text{ g CO}_2/\text{H}_2\text{O} \left(\frac{1 \text{ g NaHCO}_3}{.369 \text{ g CO}_2/\text{H}_2\text{O}} \right) = 1.19 \text{ g NaHCO}_3 \sim 1.2 \text{ g NaHCO}_3$$

$$\frac{1.19 \text{ g NaHCO}_3}{1.35 \text{ g Total}} \times 100\% = 88.1\% \sim 88\% \text{ NaHCO}_3$$

12% KCl

Trial 2

$$.54 \text{ g CO}_2/\text{H}_2\text{O} \left(\frac{1 \text{ g NaHCO}_3}{.369 \text{ g CO}_2/\text{H}_2\text{O}} \right) = .921 \text{ g NaHCO}_3$$

$$\frac{.921 \text{ g NaHCO}_3}{1.44 \text{ g Total}} \times 100\% = 63.96\% \text{ NaHCO}_3 \sim 64\% \text{ NaHCO}_3$$

36% KCl

Claim:

The mixture of $\text{KCl}/\text{NaHCO}_3$, tested at various ^{sample} portions in the mixing container yielded the following:

Trial 1 showed 88% NaHCO_3 and 12% KCl

Trial 2 showed 64% NaHCO_3 and 36% KCl

Evidence and Analysis:

From the two trials shown in the data, it is apparent that the mixture of NaHCO_3 and KCl was not mixed thoroughly. The process of heating to constant mass was followed in both experimental trials.

It is possible that results could involve error if the HCl used to clean the crucible was not entirely removed prior to testing. Any residual HCl will react with the NaHCO_3 to generate sodium chloride, water and CO_2 .

The process used in both trials ensured thorough cleaning however. It is not likely that appreciable amounts of contaminant were present.

In the first trial at one point the NaHCO_3 started to melt. Molten NaHCO_3 can potentially react with silica present in the crucible, releasing additional CO_2 and slightly affecting the outcome of the analysis. Heat was reduced immediately upon observation of molten NaHCO_3 .

Limitations of precision in drying involved using a scale ~~the~~ with capacity limited to .01 grams. Measurements were then limited to two significant figures.

Reading and Reflection:

Post-Lab questions:

If the laboratory technician accidentally mixed NaHCO_3 with $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$ in stead of KCl , results would be affected. The hydrate $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$ will release 2 moles of water for every mole of BaCl_2 present. The additional release of H_2O will result in additional mass loss when heating. As % NaHCO_3 is calculated from lost $\text{CO}_2 + \text{H}_2\text{O}$, this additional release of H_2O will result in a % NaHCO_3 calculation that is too large.

Without knowledge of exactly how much BaCl_2 was added to the mixture, accurate measurement of the % NaHCO_3 will not be possible. It will be impossible to determine how much water mass was released from the H_2O in the BaCl_2 , and how much water mass was attributed to the NaHCO_3 .

The law of conservation of mass states that total mass of materials present after a chemical reaction is the same as total mass present at the ~~end~~ beginning of the reaction. (Brown-Lemay, p. 38, 2009). The law of constant composition states that the composition, elemental, of a pure compound is always the same, no matter ~~what~~ where the sample was generated or found (Brown-Lemay, p. 8, 2009).

According to the National Center for Biotechnology Information, sodium bicarbonate (NaHCO_3) will begin decomposition at 50°C , (Center for Biotech Info, 2014), ~~NaHCO₃~~ NaHCO_3 will melt at 270°C (chemicalbook.com, 2008). Sodium carbonate, will melt at 854°C (chemicalbook.com, 2008).

The application of the law of definite composition allows for the correct analysis of NaHCO_3 from lost CO_2 and H_2O . In the lab, maintaining a temperature between 50°C and 270°C will insure NaHCO_3 decomposition without melting. The process of heating to constant mass ensures that all CO_2 and H_2O have been removed.